**Reveal® distal attachment cap**

Bouchon de fixation distale Reveal®  
Reveal® distale Befestigungskappe  
Tappo per attacco distale Reveal®  
Tapón de conexión distal Reveal®  
Tampa de fixação distal Reveal®  
Distal fastgørelseshætte fra Reveal®  
Reveal® distale bevestigingsdop  
Reveal® distal bağlantı başlığı  
Reveal® 말단 장착 캡  
Колпачок дистальной насадки Reveal®

**Reorder No. 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775**

Référence de commande 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
Nachbestell-Nr. 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
N. di riordino 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
N.º de pedido 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
N.º de encomenda 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
Genbestillingsnr. 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
Nabestelnr. 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
Yeni Sipariş Numarası: 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
주문 번호 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775  
№ повторного заказа: 00711770 / 00711771 / 00711772 / 00711773 / 00711774 / 00711775

**INSTRUCTIONS FOR USE**

MODE D'EMPLOI  
GEBRAUCHSANLEITUNG  
ISTRUZIONI PER L'USO  
INSTRUCCIONES DE USO  
INSTRUÇÕES DE UTILIZAÇÃO  
BRUGSANVISNING  
GEBRUIKSAANWIJZING  
KULLANIM TALİMATLARI  
사용 설명서  
ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ

US Endoscopy  
5976 Heisley Road  
Mentor, OH 44060  
800-548-4873  
www.steris.com  
A subsidiary of STERIS Corporation

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**This product is not made with natural rubber latex**

**Intended Use:**

The Reveal® distal attachment cap has been designed to be attached to the distal end of the endoscope to facilitate endoscopic therapy. The Reveal® distal attachment cap is intended for the following:

- Gastrointestinal mucosal resection (endoscopic mucosal resection)
- Keeping the suitable depth of endoscope's view field

Product Description	Product Number	Outer Diameter	Fits scopes
Reveal® distal attachment cap	00711770	11.35mm	8.9 - 9.5mm
Reveal® distal attachment cap	00711771	11.8mm	9.5 – 9.9 mm
Reveal® distal attachment cap	00711772	12.4mm	9.8 - 10.5mm
Reveal® distal attachment cap	00711773	13.4mm	10.9-11.8mm
Reveal® distal attachment cap	00711774	15.0mm	12.6 – 13.2mm
Reveal® distal attachment cap	00711775	15.7mm	13.4 – 13.9mm

The above information should be used as a reference guide. Endoscope specifications may vary and designs are subject to change. Please notify your local Product Specialist for questions regarding endoscope compatibility and sizing.

**Warnings and Precautions:**


- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- Do not articulate the endoscope abruptly; this could cause patient harm such as mucosal membrane injury.
- When applying the Reveal® distal attachment cap to the distal end of the endoscope, hold the endoscope as close to the distal end as possible. This will minimize the force applied to the bending section of the endoscope and prevent damage to the bending portion of the endoscope.

**Contraindications:**

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.

**Directions for Use:**

1. Prior to clinical use, inspect and familiarize yourself with the device. If damage (e.g. tears, misshapen) is noted, do not use this product and contact your local Product Specialist.
2. Read the "Instructions for Use" and review the figure below.
3. **Verify that the outer diameter of the endoscope is compatible with the inner diameter of the Reveal® distal attachment cap.** (Please refer to compatibility chart above)
4. Dry the tip of the endoscope by wiping with a 4x4 gauze, or dry wash cloth. It is important that the Reveal® distal attachment cap and tip of the endoscope are dry to ensure that the device stays attached to the endoscope. **Do not use Vaseline (or any lubricant that contains petroleum jelly), cooking oils, silicone spray alcohol, or xylocaine spray to lubricate the instrument. Doing so could cause equipment damage or cause the instrument to fall off the endoscope inside the patient.**
5. Align the side hole of the Reveal® distal attachment cap with the Objective Lens. (See figure below)

6. Place the Reveal<sup>®</sup> distal attachment cap onto the distal end of the endoscope.
7. Slide the distal end of the endoscope to the Alignment Line.
8. Secure the Reveal<sup>®</sup> distal attachment cap to the endoscope using medical grade tape ensuring **not to cover the Side Hole**.
9. Gently pull on the Reveal<sup>®</sup> distal attachment cap to ensure that it is secure.
10. Insert the endoscope into patient and use as indicated.
11. Suction the fluid that accumulates on the Objective Lens by positioning the Side Hole in contact with the mucosal membrane. This facilitates suctioning of the fluid.
12. At the end of the procedure remove the endoscope from the patient.
13. Remove the medical grade tape from the endoscope. Remove the Reveal<sup>®</sup> distal attachment cap from the endoscope and discard before reprocessing the endoscope

#### Product Disposal:

 After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Issued Date: October 2019**

#### Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit ® oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met © of ™ gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.
별도로 명시되어 있지 않는 한, ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 STERIS Corporation이 소유하고 있는 상표입니다.