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**Modelo: 00712811**

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FABRICANTE RECOMENDA USO ÚNICO.

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**FABRICANTE:**

**United States Endoscopy Group, Inc. (US Endoscopy)**

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---

*listening... and delivering solutions\**

5976 Heisley Road  
Mentor, OH 44060 USA  
**phone** +1 440 / 639.4494  
**fax** +1 440 / 639.4495  
**email** global@usendoscopy.com  
[www.usendoscopy.com](http://www.usendoscopy.com)

**This product has been manufactured not to include latex.**

**Intended Use:**

The bite block is used to protect the endoscope insertion tube from being bitten by the patient. Bite blocks are available with and without retaining straps.

**Warnings and Precautions:**

- When placing the bite block in the patient's mouth, care should be exercised to avoid placing fingers between the bite block and the patient's teeth.
- Flexible endoscopy procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to the performance of an endoscopic procedure.
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- The bite block is a disposable medical device and not intended for reuse. Any institution, practitioner, or third party, who reprocesses, refurbishes, remanufactures, resterilizes, and/or reuses this disposable device must bear full responsibility for its safety and effectiveness.
- Assess the patient's oral capacity to ensure it is compatible with the size requirements of the bite block.
- Assess patient for appropriate or intact dentition as well as oral and maxillofacial trauma prior to using and after removal of the bite block.

**Contraindications:**

Contraindications include those specific to any endoscopic procedure, as well as those specific to any endoscopic procedure requiring the use of a bite block.

**Directions for use:**

1. Prior to clinical use, familiarize yourself with the device and read all the Instructions for Use.
2. Ensure the prongs of the nasal cannula fit into the orifices in the face plate and that the nasal cannula allows passage of the endoscope through the bite block.
3. Inspect bite block and should there be evidence of any defects or damage **DO NOT USE THIS PRODUCT AND CONTACT YOUR LOCAL PRODUCT SPECIALIST.**
4. Secure retaining strap to first cleat by inserting the retaining strap through posterior aspect of side port and threading it overcleat.
5. Place the bite block in the patient's mouth.
6. Extend retaining strap around patient's head, and extend strap so as to fit securely, but comfortably.
7. Secure retaining strap by inserting the free end through posterior aspect of side port and thread over the second cleat.
8. **If O<sub>2</sub> delivery through bite block is warranted:**
  - Remove nasal cannula from patient's nares and place each prong into respective hole in bite block.
  - Adjust oxygen as needed.
  - Continue to monitor patient as would normally be done during endoscopic procedure.

**Product Disposal:**

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Issued Date: July 2013**

**Warning:**

An issued or revision date for these instructions is included for the user's information. In the event two years has elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, as US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) eingetragen wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, il design dello slogan listening...and delivering solutions di US Endoscopy e tutti i marchi contrassegnati da ® o ™ sono stati registrati presso l'US Patent and Trademark Office degli Stati Uniti o sono marchi di fabbrica di proprietà di US Endoscopy.
Listening...and delivering solutions® o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 USEndoscopy가 소유하고 있는 상표입니다.