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Rx Only (U.S.A.)

Intended Use:

- The disposable AdvanCE® – capsule endoscope delivery device is a 2.5 mm single sheathed device indicated for transendoscopic delivery of the Given PillCam SB Video Capsule to the stomach or duodenum in patients who are either unable to swallow the video capsule, or unable to pass the video capsule beyond the pylorus in sufficient time to complete the desired diagnostic evaluation.

Warnings and Precautions:

- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to complications, hazards and techniques prior to the performance of any endoscopic procedure.
- Transendoscopic capsule placement requires skill and experience in endoscopic esophageal intubation **with an accessory device seated at the distal tip of the endoscope**. The use of this device is not recommended if the clinician lacks the required experience and proficiency.
- Baseline EGD is recommended prior to transendoscopic capsule placement due to reduced luminal field of view when intubating with this device. The stomach should be decompressed upon completion of the baseline EGD to minimize gastric looping on repeat endoscopy for capsule delivery.
- Confirm engagement between the capsule cup and the stainless steel distal thread connection prior to patient intubation to minimize the risk of accidental or premature disengagement and potential patient complications.
- During intubation and prior to placement of the capsule, the delivery device's finger ring handle should be maintained in a retracted position to avoid inadvertent or premature deployment of the capsule.
- Ensure that the endoscope accessory channel is compatible with this device (2.8mm or greater).
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- This disposable medical device is not intended for reuse. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, resterilizes, and/or reuses this disposable medical device must bear full responsibility for its safety and effectiveness.

Contraindications:

The Given PillCam SB Video Capsule is indicated for patients 10 years and older. Certain patient conditions, such as surgically altered or atypical anatomy, may hinder or prevent the endoscopic placement of the capsule. Contraindications include those specific to any endoscopic procedure, as well as those contraindicated for capsule endoscopy.

Pre-Use Instructions:

Prior to clinical use you should familiarize yourself with the device.

1. Read the "Instructions for Use" and review the photographs.
2. Inspect the package for shipping or handling damage. If damage is evident, do not use this device.
3. The capsule cup is packaged separately in a small bag located along the length of the device's catheter.

Instructions for Use:

1. Prepare patient for transendoscopic capsule placement after performance of baseline EGD.
2. Open the larger poly bag, remove the paper tape restraint from the coiled device, and gently uncoil the catheter assembly.
3. Read and then remove the cardboard tag and the capsule cup bag from the distal end of the catheter by sliding it over the yellow tip protector.
4. **DO NOT REMOVE THE YELLOW TIP PROTECTOR FROM THE CATHETER. THE TIP PROTECTOR SHOULD REMAIN ON THE DISTAL END OF THE CATHETER FOR PASSAGE THROUGH THE SCOPE'S ACCESSORY CHANNEL.**
5. Remove capsule cup from smaller bag.
6. Prior to patient intubation, open the lid of the biopsy valve and pass the catheter (leaving the yellow tip protector attached) into the endoscope's accessory channel using short strokes (see Fig. 1).
7. Once the catheter assembly exits the distal end of the endoscope's accessory channel, the yellow tip protector should be removed by either tugging gently or *unthreading* in a counter-clockwise direction. The yellow tip protector should be set aside for subsequent placement back onto the stainless steel distal connection at the completion of the procedure (as described in step 18).
8. Holding the finger ring handle in a retracted position (see Fig. 2) to ensure that the deployment cable is restrained within the catheter, the capsule cup can now be threaded onto the stainless steel distal connection by gently holding the capsule cup between thumb and forefinger and *threading* in a clockwise direction (see Fig. 3) until threading stops and the catheter begins to turn. Gently tug on the capsule cup to confirm proper connection.
9. Practice advancing and retracting the finger ring handle while observing the deployment cable as it advances into and out of the capsule cup. This brief practice will further familiarize you with the device mechanics. **Note: The finger rings do not advance the full length of the handle shaft.**

10. Retract the cable into the catheter and out of the capsule cup by holding the device's finger ring handle in a retracted position (see Fig. 2). **Load the capsule into the capsule cup with the optical dome facing out (see Fig. 4), pushing the capsule in until you hear an audible click.**
11. Gently pull the proximal end of the catheter until the capsule cup is seated against the endoscope's distal tip (see Fig. 5).
12. Lubricate the outer surface of the capsule cup and endoscope in standard fashion.
13. **Caution: The finger ring handle must remain in the fully retracted position until desired deployment location is reached to avoid inadvertent or premature deployment of the capsule** (see Fig. 6).
14. After esophageal intubation is achieved, **while maintaining the finger ring handle in a fully retracted position**, the capsule cup may be advanced from the distal tip of the endoscope under direct endoscopic visualization by advancing the catheter from its proximal end. This may enhance the luminal field of view.
15. Advance the capsule cup and endoscope into the stomach, and if necessary and safely achievable, beyond the pylorus. If post-pyloric position is not achievable, deploy capsule in the stomach in the manner described in #16 below. If post-pyloric position is achieved, the duodenum should be safely insufflated and remain insufflated for capsule deployment.
16. Once desired delivery site is determined by the clinician, **deployment is achieved by deliberately advancing the finger ring handle forward until it stops** (see Fig. 7). This will advance the deployment cable through the capsule cup. **Visual confirmation of capsule delivery must be confirmed endoscopically. This is achieved by seating the capsule cup against the endoscope's lens while looking through the clear capsule cup base until the radial slits in the distal end of the capsule cup are seen. The clinician must visually appreciate the advanced cable, the empty capsule cup, and the radial slits to confirm delivery. If the capsule does not fully deploy**, retract the finger ring handle, slightly alter the position/angulation of the endoscope and repeat the process of deliberately advancing the finger ring handle forward until it stops. **Repeat visual confirmation of capsule delivery as described above.**
17. Upon completion of capsule delivery, retract the finger ring handle until the deployment cable is retracted into the catheter and out of the capsule cup. While maintaining this position, and with the capsule cup seated against the endoscope's distal tip, the endoscope may be extubated from the patient.
18. To remove the device from the endoscope, hold the finger ring handle in a retracted position to ensure that the deployment cable is restrained within the catheter and out of the capsule cup, remove the capsule cup from the distal end of the catheter by *unthreading* in a counter-clockwise direction, slide the yellow tip protector back onto the stainless steel distal connection, and withdraw the device from the endoscope's accessory channel.

Product Disposal: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: April 2010

Warning: An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

PillCam™ SB is a trademark of Given Imaging Ltd.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit © oder ™ bezeichneten Marken sind eingetragene Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) registriert wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med © eller ™ er registreret i U.S. Patent and Trademark Office eller er varemærker, der tilhører US Endoscopy
Listening...and delivering solutions®, het ontwerp *listening...and delivering solutions* van US Endoscopy en alle aanduidingen met © of ™ zijn geregistreerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청 (U.S. Patent and Trademark Office) 에 등록되어 있거나 US Endoscopy 가 소유하고 있는 상표입니다.