Exacto® cold snare
Anse à froid Exacto®
Exacto® kalte Schlinge
Ansa a freddo Exacto®
Lazo en frío Exacto®
Ansa fria Exacto®
Exacto®-koldslynge
Exacto® niet-geactiveerde poliepsnoerder
Exacto® soğuk snare
Exacto® 콜드 스네어
Петля для холодной экстракции Exacto®

Reorder No. 00711115 / 00711094 / 00711095
Référence de commande no 00711115 / 00711094 / 00711095
Nachbestell-Nr. 00711115 / 00711094 / 00711095
N. di riordino 00711115 / 00711094 / 00711095
N.º de pedido 00711115 / 00711094 / 00711095
N.º de encomenda 00711115 / 00711094 / 00711095
Genbestillingsnummer 00711115 / 00711094 / 00711095
Nabestelnr. 00711115 / 00711094 / 00711095
Yeni Sipariş Numarası. 00711115 / 00711094 / 00711095
주문 번호 00711115 / 00711094 / 00711095
№ для заказа 00711115 / 00711094 / 00711095

INSTRUCTIONS FOR USE
MODE D'EMPLOI
GEBRAUCHSANLEITUNG
ISTRUZIONI PER L'USO
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
BRUGSANVISNING
GEBRUIKSAANWIJZING
KULLANIM TALİMATLARI
사용 설명서
ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ

This product is not made with natural rubber latex.

Intended Use:

The Exacto® cold snare is intended to be used without diathermic energy for the endoscopic resection of polyp tissue in the gastrointestinal tract.

<u>Description</u>	<u>Product Number</u>	<u>Box Quantity</u>	<u>Sheath Diameter</u>	<u>Sheath Length</u>	<u>Wire Diameter</u>	<u>Approximate Snare Size (mm) Width and Length</u>
Exacto® cold snare	00711115 00711094 00711095	10 per box 20 per box 50 per box	2.4 mm	230 cm	0.30 mm (0.012 inches)	9 X 19 mm

Warnings and Precautions:

- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure.
- A thorough understanding of the technical principles, clinical applications, and risks associated with cold snaring (non-electrical) polypectomy and tissue resection is necessary before using this product.
- These devices are compatible with an endoscope channel of 2.8 mm or larger.
- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- Disengage the snare from the polyp if there is a risk of complication.
- Care should be exercised when grasping tissue to avoid inadvertently grasping tissue or organs not intended for retrieval.
- If resistance to insertion is encountered, reduce the angulation (or lower the forceps elevator if applicable) until the instrument passes smoothly.
- Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.
- The following conditions may not allow the device to function properly or cause patient injury:
 - 1) Attempting to advance the handle to the open position with too much speed or force,
 - 2) Attempting to pass or open the device in an extremely articulated endoscope,
 - 3) Attempting to actuate the device in an extremely coiled position and/or,
 - 4) Attempting to actuate the device when the handle is at an acute angle in relation to the sheath.

Contraindications:

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.

Prior to Use:

1. Prior to clinical use, inspect and familiarize yourself with the device. If there is evidence of damage, **do not use this product and contact your local Product Specialist.**
2. Ensure that the device is compatible with an endoscope channel of 2.8mm or larger before insertion.
3. Remove the device from the package and uncoil the entire device and drape in a "U" shaped configuration, holding the proximal end in one hand and the distal sheath in the opposite hand.
4. Move the finger rings back and forth to confirm that the snare loop opens and closes smoothly prior to inserting into the endoscope.
5. Ensure the snare loop is fully retracted into the catheter prior to insertion into the endoscope.

- Note: The Exacto® cold snare is *not* intended to be used with diathermic energy and, therefore, has no diathermic handle connection.

Directions for use:

1. Insert the device into the biopsy port of the endoscope using short strokes, approximately 1"-1.5" (2.5cm – 3.8cm) in length throughout device passage to avoid sheath kinking.
2. When the polyp has been endoscopically visualized, extend the distal tip of the catheter into the endoscopic field of view.
3. Gently deploy the snare loop to a fully open position.
4. Place the snare around the polyp tissue being resected using proper endoscopic technique.
5. After the polyp(s) has been satisfactorily resected, the snare should be retracted back into the catheter.
6. The polyp tissue should be collected and the specimen(s) prepared according to standard technique for histologic evaluation.

**Product Disposal:**

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: September 2019

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit © oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met ® of ™ gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.
별도로 명시되어 있지 않는 한, ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 STERIS Corporation이 소유하고 있는 상표입니다.