## **eTrap® – polyp trap**

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**eTrap®** – Polypenfalle

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800-548-4873

[www.steris.com](http://www.steris.com)

A subsidiary of STERIS Corporation

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**This product is not made with natural rubber latex.**

**Intended Use:**

The single use **eTrap® – polyp trap** is used for suction retrieval of endoscopically removed polyps.

**Warnings and Precautions:**

- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- Two specimen retrieval strainers accompany the eTrap®. Always insert the extra strainer into the eTrap® upon removal of the first strainer; this practice will guard against inadvertent suctioning of a polyp through the eTrap® and into the main suction system.

**Directions for Use:**

1. Prior to clinical use, open the package and remove the device.
2. Familiarize yourself with the device(s) and read all the Instructions for Use.
3. Inspect the package and device for shipping or handling damage, i.e. cracked clear vessel, bent, misshapen, or missing specimen retrieval strainer(s), cracked, torn, or missing connection tube. If damage is evident **DO NOT USE THESE DEVICES, SAVE THEM FOR RETURN, AND CONTACT YOUR LOCAL PRODUCT SPECIALIST.**
4. Remove both specimen retrieval strainers from the packaging and insert one of the strainers into the rectangular side port of the eTrap®.
5. Keep the second specimen retrieval strainer available for use as described in steps 10 and 11.
6. Connect the tube from the top of the eTrap® to the endoscope's suction barb.
7. Connect the main suction canister's tubing to the eTrap® device's outgoing suction barb.  
**Caution:** Hang trap ensuring that the suction canister's tubing does not kink. Kinked tubing will prevent the eTrap from functioning properly, or may decrease or eliminate the ability to suction effectively.
8. Suction through the endoscope using standard endoscopic technique until the excised polyp is entrapped in the specimen retrieval strainer.  
**Note:** The magnifying window on top of the eTrap® aids in visualization.
9. Remove the specimen retrieval strainer from the eTrap®.
10. Insert the extra specimen retrieval strainer into the eTrap® for additional usage.
11. The process of removing and replacing the specimen retrieval strainers in tandem (as needed), prevents loss of suction and allows the availability to retrieve additional polyps while the previously retrieved polyp is prepared by clinical staff.  
**Note:** Preparation of the retrieved polyp is achieved by removing it from the specimen retrieval strainer and placing it in the appropriate fixative for histopathology in accordance with institutional guidelines.
12. When the eTrap® is no longer needed, disconnect the main suction canister's tubing from the eTrap®, disconnect the eTrap® device's tubing from the endoscope and reconnect the main suction canister's tubing to the endoscope.

**Product Disposal:**

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Packaging is constructed of polyethylene terephthalate. Please recycle per your state or local recycling laws.

**Issued Date: October 2019**

**Warning:**

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit ® oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met ® of ™ gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.
별도로 명시되어 있지 않는 한, ® 또는 ™으로 표시된 모든 상표는 미국 특허청(**U.S. Patent and Trademark Office**)에 등록되어 있거나 **STERIS Corporation** 이 소유하고 있는 상표입니다.