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Modelo: 00711175

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FABRICANTE RECOMENDA USO ÚNICO.

FABRICANTE:

United States Endoscopy Group, Inc. (US Endoscopy)

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email global@usendoscopy.com

www.usendoscopy.com

This device is not made with natural rubber latex.

Intended Use:

The Talon® grasping device is intended for retrieval of foreign bodies, tissue specimens, stones or calculi in endoscopic procedures of the gastrointestinal tract.

Warnings and Precautions:

- Do not use this device if there is any apparent shipping or handling damage.
- The endoscopic retrieval of foreign objects or food bolus should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Only use this device with diagnostic and therapeutic scopes.
- Use in conjunction with an overtube is recommended.
- Use of an overtube should not be attempted unless proficiency in associated techniques has been developed by the clinician.
- Failure to use an overtube during foreign body or food bolus retrieval may lead to complications. Such complications should be monitored for during and after the procedure and treated appropriately.
- Use of overtubes has been associated with esophageal or pharyngeal laceration or perforation. Such complications should be monitored for during and after the procedure and treated appropriately.
- Avoid blindly passing the Talon® grasping device past any foreign object, particularly if the entire lumen is blocked.
- Care should be exercised when grasping the targeted object/material to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Keep traction on the device during retrieval so that the retrieved object does not become loosened, dislodged or aspirated into the trachea.
- Optimal technique to remove a foreign object or food bolus will vary depending upon the type of object/bolus or clinical conditions faced during retrieval (eg. bolus position, bolus consistency, anatomical, or disease issues such as strictures, bolus-induced ulcerations, etc.).
- The use of the Talon® grasping device or any technique described in any document or publication is at the discretion of the medical personnel attempting the retrieval.
- US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- The Talon® grasping device is intended for single patient use only. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, resterilizes, and/or reuses this disposable device must bear full responsibility for its safety and effectiveness.

Contraindications:

- Contraindications include those specific to any endoscopic procedure, endoscopic procedures to remove a foreign body or food bolus, use of an overtube, or any endoscopic procedure performed in conjunction with an overtube.
- Contraindications specific to overtube usage are: esophageal bleeding and/or laceration; laryngeal perforation; trauma to teeth, gums and/or pharynx; aspiration pneumonia.

Directions for Use:

Prior to clinical use you should familiarize yourself with the device.

1. Read the "Warnings and Precautions" and the "Directions for Use."
2. Open the package, remove the protective tip protector from the device's distal end and uncoil the entire device prior to manipulation of the handle.
3. Remove the red warning label from the handle.
4. Drape the device in a "U" shaped configuration, holding the proximal device end in one hand, and the distal sheath end in the opposite hand.
5. With the device in this position, carefully retract the handle until the grasping claws are completely retracted into the sheath.
6. Check that the accessory channel diameter of the endoscope is compatible (2.8mm or larger) with the diameter of the Talon® grasping device (2.5mm sheath).
7. If used, consult the overtube manufacturer's instructions for use separately.
8. Once the object has been endoscopically identified, advance the retracted Talon® grasping device into the accessory channel of the endoscope using short strokes until the distal end of the sheath is endoscopically visualized.
9. Open the Talon® grasping device by advancing the handle forward until it stops. Confirm that the device is fully open via endoscopic observation.
The following conditions may cause the device to function improperly: (1) attempting to pass or open the device in an extremely articulated endoscope, (2) attempting to actuate the device in an extremely coiled position and/or (3) actuating the device when the handle is at an acute angle in relation to the sheath.
10. Endoscopically manipulate the Talon® grasping device over the object or bolus.
11. Gently advance the catheter forward while retracting the handle to capture the foreign body. Continuous traction should be applied on the handle to keep the Talon® grasping device closed.
12. Keeping the distal device end endoscopically visible (as not to lose sight of the object or bolus during removal), extubate the endoscope and Talon® grasping device from the patient. The retrieved object or bolus can be removed from the device by advancing the handle forward to open the device . If used, extubate the overtube once the endoscope and Talon® grasping device have been removed from the patient.
13. Once the endoscopic retrieval has been completed, remove the device from the endoscope.

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: September 2013

Note:
An issued or revision date for these instructions is included for the user's information. In the event that more instructions for use are desired, please go to www.usendoscopy.com or contact US Endoscopy directly.

Packaging is constructed of polyethylene terephthalate. Please recycle.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind Marken, die beim U.S. Patent and Trademark Office (US-Patent- und Markenamt) eingetragen wurden, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, l'ideazione dello slogan listening...and delivering solutions di US Endoscopy e tutti i marchi contrassegnati da ® o ™ sono stati registrati presso l'US Patent and Trademark Office degli Stati Uniti o sono marchi di fabbrica di proprietà di US Endoscopy.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, het ontwerp listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.