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Modelo: 00711178

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FABRICANTE:

United States Endoscopy Group, Inc. (US Endoscopy)

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This device is not made with natural rubber latex.

Intended Use: The disposable Raptor® grasping device-mini is intended to be used to grasp tissue and/or retrieve foreign bodies, excised tissue and stents during endoscopic procedures.

Description	Product Number	Sheath Diameter	Working Length
Raptor® grasping device-mini	00711178	1.9mm	200cm

Warnings and Precautions:

- The endoscopic retrieval of foreign bodies, food bolus, or resected tissue specimens should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Use in conjunction with an overtube is recommended. Failure to use an overtube during foreign body or tissue retrieval may lead to complications.
- Care should be exercised when grasping the targeted object to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Optimal technique to remove foreign bodies will vary depending upon the type or clinical conditions faced during retrieval.
- When using the Raptor® grasping device - mini, always wear appropriate personal protective equipment. Scattering of blood, mucous or other patient debris could occur if the device is withdrawn too quickly.
- This device has not been designed nor is it intended to be used with electro-surgical devices. If use of an electro-surgical device is desired, completely remove the Raptor® grasping device - mini from the endoscope prior to actuating the electro-surgical device.
- This device was not designed to remove permanent stents.
- U.S. Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- This disposable device is not intended for reuse. Any institution, practitioner, or third-party who reprocesses, refurbishes, remanufactures, re-sterilizes, and/or reuses this disposable device must bear full responsibility for its safety and effectiveness.

Contraindications:



- Contraindications include: those specific to any endoscopic procedure, endoscopic procedures to remove a foreign body or food bolus, as well as use of an overtube, or any endoscopic procedure performed in conjunction with an overtube.

Prior to use:

1. Read the entire Instructions for Use and familiarize yourself with the device.
2. Open the package and remove the device. Before manipulating the Slider on the Handle, remove the tip protector from the device's distal end and uncoil the entire device.
3. Visually inspect the device for damage.
4. Drape the catheter in a "U" shaped configuration.
5. Actuate the device by moving the Slider back and forth to confirm that the grasping jaws open and close smoothly. (See Figure 1). *If the unit does not function properly, or there is evidence of damage (e.g. bends, kinks, or misaligned jaws) **do not use this product and contact your local Product Specialist.***

Directions for Use:

1. Identify the object to be removed via endoscopic visualization. **Do not attempt to blindly pass the Raptor® grasping device-mini past any foreign body, particularly if the entire lumen is blocked.**
2. Before inserting the device into the endoscope, ensure that the grasping jaws are in the closed position by moving the Slider toward the Thumb Ring (see Figure 2). When inserting the instrument into the endoscope, hold the Slider firmly to prevent inadvertent opening of the grasping jaws.
3. Advance the device down the accessory channel using short strokes (1.0 to 1.5 inches) to avoid kinking the catheter.
4. If resistance is felt, reduce the angulation of the endoscope or lower the forceps elevator in the endoscope until the device passes smoothly. **Do not force the device.**
5. Advance the Raptor® grasping device - mini carefully while visualizing the object to be removed. Maintain visualization of the foreign object.
6. OPEN the grasping jaws by moving the Slider away from the Thumb Ring. (See Figure 1)
7. Place the jaws around the foreign object.
8. CLOSE the grasping jaws by moving the Slider toward the Thumb Ring, and grasp the object. (See Figure 2)
9. Secure the object in the grasping jaws.
10. Apply continuous GENTLE traction while squeezing the Slider toward the Thumb Ring. Do not use excessive force or coil the sheath outside of the endoscope. This may damage the device and render it inoperable.
11. Withdraw the object together with the endoscope from the patient. **Once the foreign object is captured, do not attempt to retract the foreign object and the device into the endoscope. This could damage the endoscope and the device.**
12. Once the endoscope and the device have been withdrawn from the patient, OPEN the grasping jaws by moving the Slider away from the Thumb Ring (See Figure 1), to release the object from the grasping jaws.
13. If removal of additional foreign objects is desired, OPEN and CLOSE the device to confirm that the grasping jaws continue to operate smoothly and then repeat steps 1-13.

 **Product Disposal:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: January 2014

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ gekennzeichneten Marken sind Marken, die beim U.S. Patent and Trademark Office (US-Patent- und -Markenamt) eingetragen sind, oder Marken, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, il design dello slogan *listening...and delivering solutions* di US Endoscopy e tutti i marchi contrassegnati da ® o ™ sono stati registrati presso l'US Patent and Trademark Office degli Stati Uniti o sono marchi di fabbrica di proprietà di US Endoscopy.
Listening...and delivering solutions[®], o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com [®] ou [™] estão registados no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, designet US Endoscopy listening...and delivering solutions og alle varemærker markeret med ® eller ™ er registreret i USA's patent- og varemærkeregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™ 으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 USEndoscopy 가 소유하고 있는 상표입니다.