**Infinity® cytology brush**  
**Brosse cytologique Infinity®**  
**Infinity® Zytologiebürste**  
**Spazzolino per citologia Infinity®**  
**Cepillo de citología Infinity®**  
**Escova de citologia Infinity®**  
**Infinity® cytologibørste**  
**Infinity® cytologieborstel**  
**Infinity® sitoloji fırçası**  
**Infinity® 세포검사 브러시**  
**Цитологическая щетка Infinity®**

REORDER NO. 00711499  
RÉF. DE COMMANDE 00711499  
NACHBESTELLNR. 00711499  
N. DI RIORDINO 00711499  
Nº DE PEDIDO 00711499  
N.º DE ENCOMENDA 00711499  
GENBESTILLINGSNR. 00711499  
NABESTELNR. 00711499  
YENİ SİPARİŞ NUMARASI: 00711499  
주문 번호 00711499  
№ для ЗАКАЗА 00711499

---

**INSTRUCTIONS FOR USE**  
**MODE D'EMPOI**  
**GEBRAUCHSANLEITUNG**  
**ISTRUZIONI PER L'USO**  
**INSTRUCCIONES DE USO**  
**INSTRUÇÕES DE UTILIZAÇÃO**  
**BRUGSANVISNING**  
**GEBRUIKSAANWIJZING**  
**KULLANIM TALİMATLARI**  
**사용 설명서**  
**ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ**

---

US Endoscopy  
5976 Heisley Road  
Mentor, OH 44060  
800-548-4873  
www.steris.com  
A subsidiary of STERIS Corporation

***This device is not made with natural rubber latex.***

**Intended Use:**

The disposable **Infinity® cytology brush** is intended to be used to retrieve cytological cell samples in the gastrointestinal tract.

**Warnings and Precautions:**

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- The brush should never be forced through the endoscope's channel. If resistance is met due to excess angulation of the endoscope, it may be necessary to decrease the angulation to allow passage of the brush.

**Contraindications:**

Contraindications include those specific to any endoscopic procedure, as well as those specific to brush cytology.

**Pre-Use Instructions:**

Prior to clinical use, familiarize yourself with the device and read all the Instructions for Use. Inspect the package for shipping or handling damage. **If damage is evident, do not use this device and contact your product specialist.**

**Directions for Use:**

1. Prepare the patient for endoscopy, following recommended medical technique.
2. Remove the device from the sterile package. Uncoil the entire device and drape in a "U" shaped configuration.
3. Hold the proximal end and the distal sheath end in one hand. Remove plastic brush protector with the other hand.
4. Fully retract and deploy the spool handle and thumb ring 3-4 times to ensure the device functions properly and is not damaged. If damage is evident, do not use this device.
5. Fully retract the brush into the sheath and insert into the endoscope.
6. Advance the device into the accessory channel of the endoscope using short strokes until the distal end of the sheath is endoscopically visualized.
7. Advance the brush out of the catheter slowly using the spool handle and thumb ring and brush the desired area in order to obtain adequate cellular material.
8. Retract the brush into the sheath upon completion of brushing, in order to preserve the collection of cells.
9. Confirm the retraction of the brush into the sheath by endoscopically visualization before removing it from the endoscope.
10. Remove the brush assembly from the endoscope.
11. Properly prepare the specimen per institutional guidelines for cytological evaluation.

**Product Disposal:**

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**Issued Date:** September 2019

**Warning:** An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit ® oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met ® of ™ gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.
별도로 명시되어 있지 않는 한, ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 STERIS Corporation 이 소유하고 있는 상표입니다.