Infinity® ERCP sampling device

Dispositif d'échantillonnage Infinity® ERCP
Infinity® ERCP Probenentnahmeinstrument
Dispositivo di campionamento Infinity® ERCP
Dispositivo de muestreo Infinity® ERCP
Dispositivo de amostragem Infinity® ERCP
Infinity® ERCP-prøveudtagningsanordning
Infinity® ERCP instrument voor monsterafname
Infinity® ERCP örnek alma cihazı
Infinity® ERCP 샘플링장치
Прибор для отбора проб Infinity® ERCP

Reorder No. 00711652 / 00711653

Référence de commande 00711652 / 00711653
Nachbestell-Nr. 00711652 / 00711653
N. di riordino 00711652 / 00711653
Nº de pedido 00711652 / 00711653
N.º de encomenda 00711652 / 00711653
Genbestillingsnr. 00711652 / 00711653
Nabestelnr. 00711652 / 00711653
Yeni Sipariş Numarası 00711652 / 00711653
주문번호 00711652 / 00711653
Номер для заказа 00711652 / 00711653

INSTRUCTIONS FOR USE

MODE D'EMPLOI
GEBRAUCHSANLEITUNG
ISTRUZIONI PER L'USO
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
BRUGERVEJLEDNING
GEBRUIKSAANWIJZING
KULLANIM TALIMATLARI
사용설명서
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ

This product is not made with natural rubber latex.

Description	Product #	Sheath Diameter	Sheath Length	Approx. Brush Size (Width x Length)	Max. Guide Wire
Infinity® ERCP sampling device	00711652	3.0 mm	200 cm	4.75 mm x 29 mm	0.035in/0.89mm
Infinity® ERCP sampling device	00711653	2.5 mm	200 cm	3.0 mm x 24 mm	0.035in/0.89mm

Intended Use:

The disposable **Infinity® ERCP sampling device** is intended to be used to retrieve cytological cell samples in the gastrointestinal tract.

Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- Do not force the brush through the endoscope's channel. Reduce angulation of the scope if resistance is met.
- Do not use contrast injection with this device. Use of contrast will leak out of the short wire access hole.
- Do not rapidly remove wires from this device. Scatter of bodily fluids and debris could occur.
- When using a guide wire, hold the guidewire while brushing for cellular material. Otherwise the guidewire may move and cause patient injury.
- This device is not to be used for any purpose other than what it is intended. Use of this device for any other purpose than what it is intended could cause patient harm (eg. Perforation, cross contamination, bleeding).

Contraindications:

Contraindications include those specific to any endoscopic procedure, as well as those specific to ERCP, and any techniques associated with ductal brush cytology.

Prior to the procedure:

1. Read all the Instructions for Use and familiarize yourself with the device.
2. Open the sterile pouch and uncoil the brush.
3. Drape the catheter in a "U" shaped configuration.
4. Remove plastic brush protector.
5. Tighten the luer lock between the white handle and the gray y-port.
6. Hold the proximal device end in one hand, and the distal sheath end in the opposite hand.
7. **Rinse the brush with sterile water.**
8. Actuate the handle several times to be sure the brush functions properly. *If the unit does not function properly, or there is evidence of damage (e.g., bent brush, bent brush stem, kinked catheter), do not use this product and contact your local Product Specialist.*

Directions for Use (Long Guidewire technique):

1. Carefully retract the brush into the catheter past the radiopaque band near the tip of the catheter.
2. When using a guidewire, irrigate the guidewire lumen of the catheter with sterile water/saline.
3. After selective cannulation of the desired duct is achieved and the guidewire is placed beyond the desired site for brushing:
 - a. Advance the Infinity® ERCP sampling device over the guidewire.
 - b. With your finger, cover the hole approximately **20 cm.** from the tip of the catheter when performing long wire technique.
 - c. Observe the guide wire as it advances past the hole approximately **20 cm.** from the tip of the catheter and exits the luer lock fitting.
4. Advance the device into the accessory channel of the endoscope using short strokes (1"- 1.5") until the distal end of the sheath is endoscopically visualized.
5. Advance the brush catheter using fluoroscopic guidance to the desired location.
6. Using the spool handle and thumb ring, slowly advance the brush out of the catheter and brush the desired area in order to obtain adequate cellular material.
 - d. Brushing of the desired area will be easier with the endoscope elevator open.
7. When brushing is complete, retract the brush into the catheter and remove the Infinity® sampling device using accepted long wire exchange technique.
8. Properly prepare the specimen per institution guidelines for cytological evaluation. If salvage cytology is desired, see **Directions for Use (Salvage Cytology).**

Directions for Use (Short Guidewire technique):

1. Carefully retract the brush into the catheter past the radiopaque band near the tip of the catheter.
2. When using a guidewire, irrigate the guidewire lumen of the catheter with sterile water/saline.
3. After selective cannulation of the desired duct is achieved and the guide wire is placed beyond the desired site for brushing:
 - a. Advance the Infinity® ERCP sampling device over the guidewire.
 - b. Observe the guidewire as it exits the hole approximately **20 cm** from the tip of the catheter.
 - c. If using a guidewire locking device, **unlock** the short guidewire per manufacturer's recommendations.
 - d. Advance the device into the accessory channel of the endoscope and **relock** the guidewire.
4. Using short strokes (1" – 1.5"), continue to advance the device until the distal end of the sheath is endoscopically visualized.
5. Advance the catheter using fluoroscopic guidance to the desired location.
6. Using the spool handle and thumb ring, slowly advance the brush out of the catheter and brush the desired area in order to obtain adequate cellular material. Brushing of the desired area will be easier with the endoscope elevator open.
7. When brushing is complete, retract the brush into the catheter and remove the Infinity™ sampling device using accepted short wire exchange technique.
8. Properly prepare the specimen per institution guidelines for cytological evaluation. If salvage cytology is desired, see **Directions for Use (Salvage Cytology)**.

Directions for Use (Salvage Cytology):

1. Warning: Do not rapidly flush the catheter when performing salvage cytology. Scatter of bodily fluids and debris could occur. Cut the brush close to the bristles using an appropriate cutting device.
2. Detach brush handle from the handle luer lock.
3. Remove the entire handle assembly with the drive wire from the catheter.
4. Attach a sterile 5cc syringe filled with fresh cytology solution and attach the syringe to the handle luer lock.
5. Flush the catheter contents into the cytology container; repeat until catheter is cleared, followed by air in order to evacuate any residual contents.
6. Properly prepare the specimen per institution guidelines for cytological evaluation.

Product Disposal:

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: October 2019

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Если не указано иное, все марки, обозначенные символами © или ™, зарегистрированы в Бюро патентов и товарных знаков США, или являются товарными знаками, принадлежащими корпорации STERIS.
Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit ® oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met ® of ™ gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.