MICRO PINÇA DE BIOPSIA MORAY

Modelo: 00711889

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BRUGSANVISNING

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KULLANIM TALİMATLARI

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*listening... and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

This device is not made with natural rubber latex.

Intended Use: The disposable Moray™ micro forceps is intended to sample tissue from lesions that can occur within and outside the gastrointestinal tract (e.g. pancreas).

Description	Product Number	Sheath Diameter	Jaw Diameter	Jaw opening Width	Device Length
Moray™ micro forceps	00711889	0.80 mm	0.76 mm	4.3 mm	2300 mm

Warnings and Precautions:

- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- The sampling of tissue within and outside the gastrointestinal tract should only be performed by persons having adequate training and familiarity with endoscopic and endoscopic ultrasound techniques.
- Always use appropriate PPE when performing endoscopic and endoscopic ultrasound tissue sampling procedures.
- If using the device through an FNA needle, make sure to maintain EUS visualization of forceps at all times during usage. Losing visualization of forceps may allow the end user to advance forceps too far, potentially injuring surrounding tissue or structures.
- Care should be exercised when grasping the targeted object to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Optimal technique for tissue sampling will vary depending upon the type or clinical conditions faced during endoscopy.
- The risks associated with using this device through an FNA needle under EUS guidance are similar to the risks incurred when standard EUS FNA is performed.
- Use caution when handling included extraction pick to minimize risk of sharps injury.
- Never** swish forceps in preservative medium (i.e., formalin, methanol, alcohol fixatives, etc.) to facilitate removal of specimen from jaws if an additional pass with the forceps is warranted. Doing so may expose patient to potentially toxic material.
 - Use included extraction pick or swish in saline to facilitate specimen removal.
- If the jaws of the Moray™ micro forceps become stuck in a closed position, rinse jaws in sterile saline while actuating handle until they reopen.
 - If functionality returns to jaws, remove any excess tissue and place in specimen container.
 - If continued use is warranted continue to use the device.
 - If functionality does not return to device, do not use this product and contact your local Product Specialist for return.
 - Open up additional device if further specimen acquisition is warranted.
 - This device has not been designed nor is it intended to be used with electrosurgical devices. If use of an electrosurgical device is desired, completely remove the Moray™ micro forceps from the endoscope prior to actuating the electrosurgical device.
 - US Endoscopy did not design this device to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
 - This disposable device is not intended for reuse. Any institution, practitioner, or third party who reprocesses, refurbishes, remanufactures, resterilizes, and/or reuses this disposable device must bear full responsibility for its safety and effectiveness.

Contraindications:

- Those specific to any endoscopic procedure.
- Those specific to any endoscopic ultrasound FNA procedure or procedure where gaining access to a specific site is desired.
- Patients with known or suspected coagulopathy.

Prior to use:

- Confirm that the Moray™ micro forceps is compatible with the inner channel of the **endoscopic delivery device. (19G FNA needle, catheter, endoscope)**
- Read the entire Instructions for Use and familiarize yourself with the device.
- Place device package on flat surface with label facing up.
- Remove Tyvek® cover by grasping corner and slowly peeling it away from plastic tray. (See Figure 1)
- Remove device retention tube from tray. (See Figure 1)
- Unseat handle from tray at cutout point 1 using thumb and index finger of one hand to grasp handle. (See Figure 1)
- Maintain grasp of handle.
- Remove the rest of device by using thumb and index finger of other hand to grasp coiled sheath from cutout point 2. (See Figure 1)
- Slowly move hand to a position on the coiled sheath that is a few inches from the distal end of device. (See Figure 2)
- Uncoil the entire device slowly.
- Visually inspect the device for damage.
- Drape the sheath in a "U" shaped configuration.
- Actuate the device by moving the Slider back and forth to confirm that the jaws open and close smoothly. (See Figure 2). *If the unit does not function properly, or there is evidence of damage (e.g. bends, kinks, or misaligned jaws) do not use this product and contact your local Product Specialist.*

Directions for Use:

1. Identify the lesion to be sampled via endoscopic or endoscopic ultrasound visualization, and **maintain visualization throughout use of the Moray™ micro forceps.**
2. Place endoscopic delivery device distal end into or near target lesion, while maintaining visualization for both the lesion and delivery device distal end.
3. Keep the grasping jaws in the closed position by moving the Slider proximally towards the Thumb Ring when inserting the Moray™ micro forceps into the endoscopic delivery device.
4. Hold the Slider firmly back towards thumb ring to prevent inadvertent opening of the grasping jaws when inserting the instrument into the endoscopic delivery device.
5. Advance the Moray™ micro forceps down the endoscopic delivery device using short strokes (1.0 to 1.5 inches) to avoid kinking the sheath.
6. Reduce the angulation of the endoscope or lower the forceps elevator in the endoscope if resistance is felt. **Do not force the Moray™ micro forceps through the endoscopic delivery device.**
7. Advance the Moray™ micro forceps carefully through the endoscopic delivery device while **visualizing the lesion to be sampled.**
8. Maintain EUS or endoscopic visualization of the endoscopic delivery device within or above the lesion to be sampled.
9. Advance the Moray™ micro forceps out of the endoscopic delivery device to desired area within or above lesion.
(See Figure 3)
Caution: If using device through an FNA needle, make sure to maintain EUS visualization of forceps at all times during usage. Losing visualization of forceps may allow the end user to advance forceps too far, potentially injuring surrounding tissue or structures.
10. OPEN the jaws by moving the Slider distally from the Thumb Ring in a controlled fashion. (See Figure 4)
Caution: Forcibly moving the Slider may allow it to travel further than required causing the coiled sheath to expand, extend and then retract. This can cause device malfunction and prevent specimen capture.
11. Advance the Moray™ micro forceps with jaws open to capture the desired tissue while **maintaining EUS or endoscopic visualization.** (See Figure 5)
12. CLOSE the jaws by moving the Slider proximally toward the Thumb Ring, and slowly grasp the tissue. (See Figure 6)
Caution: Use only enough pressure to bite the tissue. Over-exertion could cause forceps to become misaligned or to fail.
13. While maintaining pressure on Slider and Thumb Ring with one hand, grasp the sheath of the Moray™ micro forceps (just above insertion point of the delivery device) and sharply, pull in a proximal direction to remove tissue sample from site.
Caution when using through 19 gauge FNA needle: The jaws of the Moray™ micro forceps may grasp a large amount of tissue allowing the forceps/sample to pull on the distal end of the needle causing it to lose position within the lesion. To help mitigate this occurrence it is critical that the clinician performing the procedure maintain control of the needle's handle during this maneuver.
14. To remove the Moray™ micro forceps from the delivery device, apply continuous GENTLE traction while squeezing the Slider proximally towards the Thumb Ring.
15. Remove the Moray™ micro forceps from endoscopic delivery device using short continuous strokes to minimize the scatter of bioburden.
16. Transfer specimen from the jaws of the Moray™ micro forceps to specimen jar and medium for preservation for pathology per institutional guidelines.
 - a. The included extraction pick (see Figure 7) may be used to assist in specimen removal if deemed necessary.
 - b. If using the extraction pick:
 - i. Unseat extraction pick from tray at cutout point 2 by grasping device with thumb and forefinger of one hand.
 - ii. Remove tip protector from extraction pick to expose tip.
 - iii. Examine extraction pick for any damage.
 - iv. Rest open jaws of forceps onto medium of choice (i.e., filter paper, non-adherent dressing, or slide)
 - v. Insert needle tip of extraction pick into the proximal aspect of the jaw cup just behind or below specimen.
 - vi. Extract specimen from cup carefully and onto medium of choice by pushing it distally.
 - vii. Repeat steps iv. through vi. if more specimens are excised.
 - viii. Dispose of extraction pick into sharps container once it is no longer needed at end of procedure.
17. You may also rinse the jaws of the Moray™ micro forceps in sterile saline while actuating handle until tissue/specimen is removed. **Caution:** Never swish forceps in preservative medium (i.e., formalin, methanol, alcohol fixatives, etc.) to facilitate removal of specimen from jaws if an additional pass with the forceps is warranted. Doing so may expose patient to potentially toxic material.
 - a. Transfer any remaining tissue/specimen to specimen container.
18. If another tissue sample is desired, visually inspect the forceps for damage and then repeat steps 3-17.



Figure 1

Thumb Ring

Slider


Figure 2

EUS Visualization:


Figure 3


Figure 4


Figure 5


Figure 6

Figure 7

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Packaging is constructed of polyethylene terephthalate. Please recycle.

Issued Date: December 2015

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das US Endoscopy „listening...and delivering solutions“-Design und alle mit ® oder ™ bezeichneten Marken sind beim U.S. Patent und Trademark Office (US-Patent- und Markenamt) eingetragene Marken oder Warenzeichen, die sich im Besitz von US Endoscopy befinden.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions design og alle varemærker markeret med ® eller ™ er registreret i USA's patentregister (U.S. Patent and Trademark Office) eller er varemærker, der tilhører US Endoscopy.
Listening...and delivering solutions®, het ontwerp Listening...and delivering solutions van US Endoscopy en alle aanduidingen met ® of ™ zijn gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 USEndoscopy가 소유하고 있는 상표입니다.