**Guardus® overtube – esophageal**

Overtube Guardus® - œsophagien  
Guardus® Overtube – Ösophagus  
Guardus® - overtube - esofageo  
Guía esofagal Guardus®  
Overtube Guardus® – esofágico  
Guardus® overtube – øsofagus  
Guardus®-overtube – oesofageaal  
Guardus® overtube – özofajeal  
Guardus® 오버튜브 - 식도  
Guardus® - тубус для обследования пищевода

**Reorder No: 00711146 / SI711146 / 00711147 / SI711147**

Référence de commande 00711146 / SI711146 / 00711147 / SI711147  
Nachbestell-Nr. 00711146 / SI0071146 / 00711147 / SI711147  
N. di riordino 00711146 / SI0071146 / 00711147 / SI711147  
N.º de pedido 00711146 / SI0071146 / 00711147 / SI711147  
N.º de encomenda 00711146 / SI0071146 / 00711147 / SI711147  
Genbestillingsnr. 00711146 / SI0071146 / 00711147 / SI711147  
Nabestelnr.00711146 / SI0071146 / 00711147 / SI711147  
Yeni Sipariş Numarası: 00711146 / SI0071146 / 00711147 / SI711147  
주문 번호 00711146 / SI0071146 / 00711147 / SI711147  
№ для заказа 00711146 / SI0071146 / 00711147 / SI711147

**Guardus® overtube – gastric**

Overtube Guardus® - gastrique  
Guardus® Overtube – Magen  
Guardus® - overtube - gastrico  
Guía gástrica Guardus®  
Overtube Guardus® – gástrico  
Guardus® overtube – gastrisk  
Guardus® overtube – voor maag  
Guardus®-overtube – gastrik  
Guardus® 오버튜브 - 위장  
Guardus® - тубус для обследования желудка

**Reorder No. 00711148 / SI711148 / 00711149 / SI711149**

Référence de commande 00711148 / SI711148 / 00711149 / SI711149  
Nachbestell-Nr. 00711148 / SI711148 / 00711149 / SI711149  
N. di riordino 00711148 / SI711148 / 00711149 / SI711149  
N.º de pedido 00711148 / SI711148 / 00711149 / SI711149  
N.º de encomenda 00711148 / SI711148 / 00711149 / SI711149  
Genbestillingsnr. 00711148 / SI711148 / 00711149 / SI711149  
Nabestelnr. 00711148 / SI711148 / 00711149 / SI711149  
Yeni Sipariş Numarası: 00711148 / SI711148 / 00711149 / SI711149  
주문 번호 00711148 / SI711148 / 00711149 / SI711149  
№ для заказа 00711148 / SI711148 / 00711149 / SI711149

---

**INSTRUCTIONS FOR USE**

Instructions d'utilisation  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Instruções de Utilização  
Brugsanvisning  
Gebruiksaanwijzing  
Kullanım talimatları  
사용 설명서  
Инструкция по применению

---

*This product is not made with natural rubber latex.*

Reorder No.	Handle Color	Device length	Inner tube ID	Outer tube ID	Outer tube OD	Endoscope OD range
00711146 / SI711146	Yellow	25cm	8.6mm	16.7mm	19.5mm	8.6mm-10.0mm
00711147 / SI711147	Orange	25cm	10.2mm	16.7mm	19.5mm	9.8mm-11.7mm
00711148 / SI711148	Yellow	50cm	8.6mm	16.7mm	19.5mm	8.6mm-10.0mm
00711149 / SI711149	Orange	50cm	10.2mm	16.7mm	19.5mm	9.8mm-11.7mm

**Overtube assembly contains (1) inner tube with tapered tip and (1) outer tube with tethered insufflation cap.**

**Intended Use:**

- The disposable overtube is a device used in conjunction with a flexible endoscope for foreign body or tissue retrieval, and/or for endoscopic procedures requiring multiple endoscope intubations.

**Warnings and Precautions:**

- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to complications, hazards and techniques prior to the performance of any endoscopic procedure.
- Use of an overtube, as well as foreign body retrieval, should not be attempted unless proficiency in technique has been developed by the clinician.
- A baseline EGD is recommended prior to overtube usage to evaluate whether there are any contraindications to using the overtube and to minimize potential patient complications.
- Placement of overtube should not be attempted with side-viewing or tangentially viewing endoscopes.
  - If use of overtube with these endoscopes is necessary, the overtube must first be placed with forward viewing endoscope (with compatible dimensions).
- Use of overtubes has been associated with esophageal or pharyngeal laceration or perforation. Such complications should be monitored for during and after the procedure and treated appropriately.
- **The Guardus® overtube is a two tube system consisting of a tapered tip inner tube and an outer tube. The tapered tip inner tube and the outer tube must be assembled together for the initial placement and repositioning of the device.**
- **Do not place or reposition the outer tube without using the tapered tip inner tube and the outer tube together. Failing to do so increases the chance of esophageal or pharyngeal perforation or laceration.**
- The true OD (outer diameter) of endoscope insertion tubes can vary from listed specifications, and therefore some endoscopes are labeled as compatible with both our standard and therapeutic overtubes.
  - **Ensure that the OD of the endoscope is compatible with the ID (inner diameter) of the inner tube with tapered tip prior to use by:**
    - **Referring to the Endoscope Reference charts with product number and scope specifications**
    - **Following the Pre-Use Instructions below**
- Do not use the device if the OD of the endoscope is incompatible with the ID of the tapered tip inner tube. **If the device is incompatible with the endoscope, consult your STERIS Endoscopy sales representative.**
- Performance may be compromised if proper lubrication is not applied to both the inner and outer surfaces of the overtube assembly components and endoscope.
- Do not insert the overtube assembly into the patient when the tapered tip section of the inner tube is in direct contact with the distal bending section of the endoscope. Performance of both the endoscope and overtube may be compromised as the distal bending section has a larger diameter than the insertion tube.
  - This could cause a noticeable bump effect that may be felt during insertion, allowing the endoscope to inadvertently advance further into the patient.
  - In order to prevent the overtube assembly from being placed when the tapered tip section on the inner tube is in direct contact with the distal bending section of the endoscope, the physician should reference where the insertion tube depth markings are in relation to the overtube markings (see Fig. 4).
- Once the overtube assembly is in place and the tapered tip inner tube has been removed, do not advance the outer tube beyond the endoscope. Any migration or advancement of the outer tube could result in tissue entrapment and mucosal injury. If repositioning is required, please follow directions as listed in steps 1-7 in the “Repositioning the Overtube” section.
- Do not attempt to reuse, reprocess, refurbish, remanufacture, or resterilize this device. STERIS Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).

**Caps/Cap-like devices:**

- The overtube assembly’s tethered insufflation cap **must be removed** to facilitate passage of caps and cap-like devices that could potentially be applied to the distal end of the endoscope. **The clinician must ensure size compatibility between specific manufacturers’ cap and cap-like device diameters and the ID (inner diameter) of the outertube (16.7 mm) prior to concomitant use.** Be aware that the overtube must be successfully placed according to these Instructions for Use prior to the use of a cap or cap-like device.

**Contraindications:**

- Esophageal bleeding and/or laceration; laryngeal perforation; trauma to teeth, gums and/or pharynx; aspiration pneumonia.
- Contraindications also include those specific to any endoscopic procedure, use of an overtube and/or foreign body retrieval.

#### Pre-Use Instructions:

Prior to clinical use, familiarize yourself with the device.

1. Read the “Instructions for Use” and review the photographs.
2. Inspect the package for shipping or handling damage (i.e. deformed tubes, holes in packaging). If damage is evident, do not use this device and **CONTACT your Product Specialist** and save the device and any accompanying materials for return.
3. Ensure that the OD (outer diameter) of the endoscope is compatible with the ID (inner diameter) of the inner tube:
  - See product number color and size specifications prior to opening package (For user convenience, the inner tube handles have been color-coded).
    - The **yellow capped** inner tube generally fits endoscope OD (outer diameter) range **8.6mm-10.0mm**.
    - The **orange capped** inner tube generally fits endoscope OD (outer diameter) range **9.8mm-11.7mm**.
  - Confirm that the tapered tip of the overtube assembly forms a minimally *gapped interface* with the OD of the endoscope, and that it slides freely over the endoscope.
4. This device is packaged in two pieces – the *tapered tip inner tube and the outer tube with tethered insufflation cap*. Prior to use, remove lubrication tag from inner tube and outer tube. **Generously lubricate both components with a water soluble lubricant** as described below in the Directions for Use. **Do not assemble without generous lubrication.**
5. **The tapered tip inner tube and the outer tube with tethered insufflation cap always must be used together during initial placement of device in the patient, or any subsequent repositioning of the device.**

#### Directions for Use:

1. Prepare patient for endoscopy. It is recommended that a minimum 60FR bite block be used to prevent inadvertent collapse of the overtube by the patient’s gums or dentition.
2. After removing the components of the overtube assembly from the packaging, **generously lubricate the inner and outer surfaces of both tubes with a water soluble lubricant.**
3. Remove the tethered insufflation cap from the outer tube. (See Fig. 1)
4. Load the **well lubricated** inner tube *into* the **well lubricated** outer tube. (See Fig. 2)
5. Lock the inner tube into the outer tube by engaging the cap threads into the handle and *twisting the cap clockwise*. (See Fig. 3)
6. Once locked, visually inspect the distal insertion end of the overtube assembly to ensure that the tapered tip forms a *gap-free interface* with the ID of the outer tube.
7. **Generously lubricate** the lower one third of the distal (insertion end) of the endoscope.
8. *Backload the well lubricated overtube assembly onto the well lubricated endoscope*, sliding it to the proximal end of the endoscope’s insertion tube so that it rests against the control head. (See Fig. 4)
9. Observe that the tapered tip of the overtube assembly forms a minimally *gapped interface* with the OD of the endoscope, and that it slides freely over the endoscope. (See Fig. 5)
10. Observe that the tapered tip of the overtube assembly is proximal to the distal bending section of the endoscope.
11. Perform a baseline EGD.
12. If the clinician determines that passage of the overtube assembly is achievable or desired, **gently raise the patient’s chin.**
13. *Gently pass* the overtube assembly *forward* over the **well lubricated** endoscope insertion tube until it reaches the desired anatomical location.
14. Observe the depth of insertion of the overtube assembly by noting the centimeter markings on the outer tube. **Avoid delivering the overtube assembly over the distal bending section of the endoscope.**
15. **Maintain the patient’s raised chin position** once the overtube assembly is in the desired location.
16. Turn the proximal lock on the color coded inner tube cap *counter-clockwise* to unlock it from the handle of the outer tube.
17. Maintain position of the outertube’s handle manually.
18. *Withdraw both the tapered tip inner tube and endoscope simultaneously* from the patient. (See Fig. 6)
19. Add additional lubrication to the lower one third of the distal (insertion end) of the endoscope and remove the inner tube from the endoscope in one smooth motion.
20. Lock the *tethered insufflation cap* onto the outer tube by twisting it *clockwise*. (See Fig. 7)
21. Once the endoscope is removed from the inner tube, the endoscopic procedure may be performed through the remaining outer tube. (See Fig. 8)

#### Repositioning the Overtube:

1. If repositioning of the outertube deeper into the gastrointestinal lumen is required, extubate the endoscope from the outer tube.
2. Remove the tethered insufflation cap.
3. **Lubricate the inner and outer surfaces of the inner tube generously.**
4. *Backload* the inner tube with tapered tip onto the **well lubricated** endoscope.
5. Extend the endoscope tip approximately 15 cm (6 inches) past the distal end of the inner tube during repeat intubation to allow for direct endoscopic visualization of the mucosal lining at the distal end of the previously placed outertube.
  - a. The clinician can then endoscopically clear the area (through insufflation) of any enveloped mucosa prior to advancing the inner tube.
6. Once the clinician is sure that the ID of the outertube is free of mucosa, the inner tube may be fully advanced and engaged with the outer tube to create a gap-free interface. The clinician can then reposition the overtube assembly to the desired position.
7. Once the desired overtube assembly position is achieved, the clinician should repeat steps 15 through 20 of the Directions for Use.

**Product Disposal:**

After use, this product may be a potential biohazard which presents a risk of cross-contamination. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Packaging is constructed of polyethylene terephthalate. Please recycle per your state or local recycling laws.

**Issued Date:** September 2019

**Warning:** An issued or revision date for these instructions is included for the user's information. In the event two years has elapsed between this date and product use, the user should contact STERIS to determine if additional information is available.

Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or are trademarks owned by STERIS Corporation.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Если не указано иное, все марки, обозначенные символами © или ™, зарегистрированы в Бюро патентов и товарных знаков США, или являются товарными знаками, принадлежащими корпорации STERIS.
Sauf indication contraire, toutes les marques accompagnées des symboles ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office, ou sont des marques commerciales détenues par STERIS Corporation.
Sofern nicht anders angegeben, sind alle mit ® oder ™ versehenen Marken beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder sind Eigentum der STERIS Corporation.
A menos que de outra forma indicado, todas as marcas com ® ou ™ são registadas na U.S. Patent and Trademark Office ou são marcas registadas propriedade da STERIS Corporation.
Medmindre andet er angivet er alle mærker markeret med ® eller ™ registreret hos U.S. Patent and Trademark Office eller er varemærker ejet af STERIS Corporation.
Tenzij anders is aangegeven, zijn alle aanduidingen met ® of ™ gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn deze aanduidingen handelsmerken van STERIS Corporation.